

Normal Hemodynamic Parameters and Laboratory Values

Normal Hemodynamic Parameters – Adult

PARAMETER	EQUATION	NORMAL RANGE
Arterial Oxygen Saturation (SaO ₂)		95 - 100%
Mixed Venous Saturation (SvO ₂)		60 - 80%
Central Venous Oxygen Saturation (ScvO ₂)		70%
Arterial Blood Pressure (BP)	Systolic (SBP) Diastolic (DBP)	100 - 140 mmHg 60 - 90 mmHg
Mean Arterial Pressure (MAP)	$SBP + (2 \times DBP)/3$	70 - 105 mmHg
Right Atrial Pressure (RAP)/ Central Venous Pressure (CVP)		2 - 6 mmHg
Right Ventricular Pressure (RVP)	Systolic (RVSP) Diastolic (RVDP)	15 - 30 mmHg 2 - 8 mmHg
Pulmonary Artery Pressure (PAP)	Systolic (PASP) Diastolic (PADP)	15 - 30 mmHg 8 - 15 mmHg
Mean Pulmonary Artery Pressure (MPAP)	$PASP + (2 \times PADP)/3$	9 - 18 mmHg
Pulmonary Artery Occlusion Pressure (PAOP)		6 - 12 mmHg
Left Atrial Pressure (LAP)		4 - 12 mmHg
Cardiac Output (CO)	$HR \times SV/1000$	4.0 - 8.0 L/min
Cardiac Index (CI)	CO/BSA	2.5 - 4.0 L/min/m ²
Stroke Volume (SV)	$CO/HR \times 1000$	60 - 100 mL/beat
Stroke Volume Index (SVI)	$CI/HR \times 1000$	33 - 47 mL/m ² /beat
Stroke Volume Variation (SVV)	$SV_{max} - SV_{min}/SV_{mean} \times 100$	10 - 15%
Systemic Vascular Resistance (SVR)	$80 \times (MAP - RAP)/CO$	800 - 1200 dynes - sec/cm ⁵
Systemic Vascular Resistance Index (SVRI)	$80 \times (MAP - RAP)/CI$	1970 - 2390 dynes - sec/cm ⁵ /m ²
Pulmonary Vascular Resistance (PVR)	$80 \times (MPAP - PAOP)/CO$	<250 dynes - sec/cm ⁵
Pulmonary Vascular Resistance Index (PVRI)	$80 \times (MPAP - PAOP)/CI$	255 - 285 dynes - sec/cm ⁵ /m ²
Left Ventricular Stroke Work (LVSW)	$SI \times MAP \times 0.0144$	8 - 10 g/m ²
Left Ventricular Stroke Work Index (LVSWI)	$SVI \times (MAP - PAOP) \times 0.0136$	50 - 62 g/m ² /beat
Right Ventricular Stroke Work (RVSW)	$SI \times MAP \times 0.0144$	51 - 61 g/m ²
Right Ventricular Stroke Work Index (RVSWI)	$SVI \times (MPAP - CVP) \times 0.0136$	5 - 10 g/m ² /beat



Edwards

Normal Hemodynamic Parameters – Adult

PARAMETER	EQUATION	NORMAL RANGE
Coronary Artery Perfusion Pressure (CPP)	Diastolic BP - PAOP	60 - 80 mmHg
Right Ventricular End-Diastolic Volume (RVEDV)	SV/EF	100 - 160 mL
Right Ventricular End-Diastolic Volume Index (RVEDI)	RVEDV/BSA	60 - 100 mL/m ²
Right Ventricular End-Systolic Volume (RVESV)	EDV - SV	50 - 100 mL
Right Ventricular Ejection Fraction (RVEF)	SV/EDV x 100	40 - 60%
Arterial Oxygen Content (CaO ₂)	$(0.0138 \times \text{Hgb} \times \text{SaO}_2) + 0.0031 \times \text{PaO}_2$	16 - 22 mL/dL
Venous Oxygen Content (CvO ₂)	$(0.0138 \times \text{Hgb} \times \text{SvO}_2) + 0.0031 \times \text{PvO}_2$	15 mL/dL
A - V Oxygen Content Difference (C(a - v)O ₂)	CaO ₂ - CvO ₂	4 - 6 mL/dL
Oxygen Delivery (DO ₂)	CaO ₂ x CO x 10	950 - 1150 mL/min
Oxygen Delivery Index (DO ₂ I)	CaO ₂ x CI x 10	500 - 600 mL/min/m ²
Oxygen Consumption (VO ₂)	C(a - v)O ₂ x CO x 10	200 - 250 mL/min
Oxygen Consumption Index (VO ₂ I)	C(a - v)O ₂ x CI x 10	120 - 160 mL/min/m ²
Oxygen Extraction Ratio (O ₂ ER)	$(\text{CaO}_2 - \text{CvO}_2)/\text{CaO}_2 \times 100$	22 - 30%
Oxygen Extraction Index (O ₂ EI)	$(\text{SaO}_2 - \text{SvO}_2)/\text{SaO}_2 \times 100$	20 - 25%

Normal Blood Laboratory Values

TEST	CONVENTIONAL UNITS (Reference Values*)	SI UNITS
Hematocrit (Hct)	Males: 42% - 52%	0.42 - 0.52
	Females: 36% - 48%	0.36 - 0.48
Hemoglobin (Hgb)	Males: 12.4 - 17.4 g/dL	124 - 174 g/L
	Females: 11.7 - 16 g/dL	117 - 160 g/L
Lactate	0.93 - 1.65 mEq/L	0.93 - 1.65 mmol/L

SI Units = International Units

**Reference Values vary by regional laboratory techniques and methods.*

Edwards is a trademark of Edwards Lifesciences Corporation. Edwards Lifesciences and the stylized E logo are trademarks of Edwards Lifesciences Corporation and are registered in the United States Patent and Trademark Office ©2009 Edwards Lifesciences LLC. All rights reserved. AR04313



Edwards Lifesciences

Edwards Lifesciences LLC

One Edwards Way · Irvine, CA 92614 USA · 949.250.2500 · 800.424.3278
www.edwards.com